

```
EEEEEEEEEE XX XX AAAAAA MM MM PPPPPPPP LL EEEEEEEEEE
EEEEEEEEEE XX XX AAAAAA MM MM PPPPPPPP LL EEEEEEEEEE
EEEEEEEEEE XX XX AAAAAA MM MM PPPPPPPP LL EEEEEEEEEE
EE XX XX AA AA MMMM MMMM PP PP LL EE
EE XX XX AA AA MMMM MMMM PP PP LL EE
EE XX XX AA AA MM MM PP PP LL EE
EE XX XX AA AA MM MM PP PP LL EE
EEEEEEEE XX XX AA AA MM MM PPPPPPPP LL EEEEEEEEE
EEEEEEEE XX XX AA AA MM MM PPPPPPPP LL EEEEEEEEE
EEEEEEEE XX XX AA AA MM MM PPPPPPPP LL EEEEEEEEE
EE XX XX AAAAAAAAAA MM MM PP LL EE
EE XX XX AAAAAAAAAA MM MM PP LL EE
EE XX XX AAAAAAAAAA MM MM PP LL EE
EE XX XX AA AA MM MM PP LL EE
EE XX XX AA AA MM MM PP LL EE
EEEEEEEEEE XX XX AA AA MM MM PP LLLLLLLLLL
EEEEEEEEEE XX XX AA AA MM MM PP LLLLLLLLLL
EEEEEEEEEE XX XX AA AA MM MM PP LLLLLLLLLL
```

```

LL      AAAAAA  BBBB BBBB  IIIIII  000000  SSSSSSSS  TTTTTTTTTT  AAAAAA  TTTTTTTTTT
LL      AAAAAA  BBBB BBBB  IIIIII  000000  SSSSSSSS  TTTTTTTTTT  AAAAAA  TTTTTTTTTT
LL      AA      AA  BB      BB  II      00      00  SS      TT      AA      AA  TT
LL      AA      AA  BB      BB  II      00      00  SS      TT      AA      AA  TT
LL      AA      AA  BB      BB  II      00      00  SS      TT      AA      AA  TT
LL      AA      AA  BB      BB  II      00      00  SS      TT      AA      AA  TT
LL      AA      AA  BBBB BBBB  II      00      00  SSSSSS  TT      AA      AA  TT
LL      AA      AA  BBBB BBBB  II      00      00  SSSSSS  TT      AA      AA  TT
LL      AAAAAAAAAA  BB      BB  II      00      00  SS      TT      AAAAAAAAAA  TT
LL      AAAAAAAAAA  BB      BB  II      00      00  SS      TT      AAAAAAAAAA  TT
LL      AA      AA  BB      BB  II      00      00  SS      TT      AA      AA  TT
LL      AA      AA  BB      BB  II      00      00  SS      TT      AA      AA  TT
LLLLLLLLLLL  AA      AA  BBBB BBBB  IIIIII  000000  SSSSSSSS  TT      AA      AA  TT
LLLLLLLLLLL  AA      AA  BBBB BBBB  IIIIII  000000  SSSSSSSS  TT      AA      AA  TT

```



```

FFFFFFFFFF  000000  RRRRRRRR
FFFFFFFFFF  000000  RRRRRRRR
FF      00      00  RR      RR
FF      00      00  RR      RR
FF      00      00  RR      RR
FF      00      00  RR      RR
FFFFFFFFFF  00      00  RRRRRRRR
FFFFFFFFFF  00      00  RRRRRRRR
FF      00      00  RR      RR
FF      00      00  RR      RR
FF      00      00  RR      RR
FF      00      00  RR      RR
FF      00      00  RR      RR
FF      000000  RR      RR
FF      000000  RR      RR

```

LB

C
C
C
4

C
C
C
5

5

C
C
C
6

C
C
C
7

C
C
C
8

.

File: LABIOSTAT.FOR
Version 'V04-000'

```
*****
*
*  COPYRIGHT (c) 1978, 1980, 1982, 1984 BY
*  DIGITAL EQUIPMENT CORPORATION, MAYNARD, MASSACHUSETTS.
*  ALL RIGHTS RESERVED.
*
*  THIS SOFTWARE IS FURNISHED UNDER A LICENSE AND MAY BE USED AND COPIED
*  ONLY IN ACCORDANCE WITH THE TERMS OF SUCH LICENSE AND WITH THE
*  INCLUSION OF THE ABOVE COPYRIGHT NOTICE. THIS SOFTWARE OR ANY OTHER
*  COPIES THEREOF MAY NOT BE PROVIDED OR OTHERWISE MADE AVAILABLE TO ANY
*  OTHER PERSON. NO TITLE TO AND OWNERSHIP OF THE SOFTWARE IS HEREBY
*  TRANSFERRED.
*
*  THE INFORMATION IN THIS SOFTWARE IS SUBJECT TO CHANGE WITHOUT NOTICE
*  AND SHOULD NOT BE CONSTRUED AS A COMMITMENT BY DIGITAL EQUIPMENT
*  CORPORATION.
*
*  DIGITAL ASSUMES NO RESPONSIBILITY FOR THE USE OR RELIABILITY OF ITS
*  SOFTWARE ON EQUIPMENT WHICH IS NOT SUPPLIED BY DIGITAL.
*
*****
```

Program LABIO STATUS

This is a utility routine for the LABIO system. It displays the status of all 16 channels of the A/D. It assumes that the terminal is a VT52 or an equivalent, e.g VT100 in VT52 mode. The display is update once every 1-9 seconds. Default is one second. There are 5 commands associated with the program

C - display status of 16 channels
P - display status by process PID
H - display help frame (timeouts after 1 min.)
E - Exit to VMS DCL
Digit(1-9) Change cycle time.

The key pad can also be used to enter commands. The special function keys on the VT52 or VT100 correspond to the first 4 commands (3 on VT52).

Typing ANY key will cause a display refresh.

Include 'LABCHNDEF.FOR'

Character*10 STATUS(4)
Character*8 XTIME
Character*9 XDATE
Parameter COMMAND_MAX = 4
Character*1 COMMAND,COMMAND_TABLE(COMMAND_MAX,2),ESCAPE,TERMINATOR
Character*63 COMMAND_DEV

External \$\$\$ NOTRAN, \$\$\$ NORMAL, \$\$\$ PARTESCAPE
External IOSM_CVTLOW, IOSM_NOECHO, IOSM_TIMED, IOS_READVBLK, IOSM_PURGE


```

Logical SUCCESS,SYSSQIOW,SYSS$ASSIGN
Integer CHANNEL,DISPLAY_FLAG,OLD_DISPLAY,COMMAND_CHAN
Integer DEF_TIME_OUT,TIME_OUT
Byte ERASE_SCREEN(2),HOME(2),ERASE_LINE(2),VT52_MODE(7)
Integer*2 IO_STATUS(4),CHAR_COUNT
Equivalence (ESCAPE,HOME),(CHAR_COUNT,IO_STATUS(2))

```

VT52 control ESCAPE Sequences

```

Data HOME,ERASE_SCREEN,ERASE_LINE
1 /'33'O,'H','33'O,'J','33'O,'K'/

```

VT100 control ESCAPE sequences

This ESC seq places a VT100 in VT52 mode

```

Data VT52_MODE/'33'O,'L','?','2','L','33'O,'J'/

```

```

Data STATUS/'Unknown','Inactive','Active'/'/
Data COMMAND_TABLE/'C','P','E','H','P','Q','S','R'/
Data DISPLAY_FLAG,ERASE_FLAG /1,TRUE./
Data DEF_TIME_OUT /1/

```

Map to the GLOBAL DATA section created by the I/O program

```

Call LABIO_INIT(0)

```

Place VT100's in VT52 mode

```

Type 500, VT52_MODE

```

Initialize Command input channel

We will read the command via a QIOW with a 1 sec timeout
 Commands are single character, to simplify matters we will
 read with no echo and convert lower to upper case.

```

Call SYSS$ASSIGN( 'TT',COMMAND_CHAN,...)
QIO_READ = %Loc(IOSM_NOECHO) + %Loc(IOSM_CVTLOW) + %Loc(IOSM_TIMED)
1 + %Loc(IOS_READVBLK)
TT_PURGE = %Loc(IOSM_PURGE)
Go To 25 ! Display Something

```

Get a command from the user, but only wait a short time (TIME_OUT)
 so we can update the screen. The input buffer is purged if a command
 was decode on the last read. (Prevents unnecessary erase loops)

```

20 DISPLAY_FLAG = OLD_DISPLAY !Default is last display
   TIME_OUT = DEF_TIME_OUT !Default time out
21 TABLE_INDEX = 1 !Assume no escape sequence
22 Call SYSS$QIOW(%Val(COMMAND_CHAN),%Val(QIO_READ+PURGE),
   1 IO_STATUS,...,%Ref(COMMAND),%Val(1),%Val(TIME_OUT),...,)
   PURGE = 0

```

```

! If escape seq., set command table pointer to second table and
! get character following escape.
  TERMINATOR = Char( IO_STATUS(3) )
  If( TERMINATOR .ne. ESCAPE ) Go To 23
  TABLE_INDEX = 2
  Go To 22      !Get char following escape

23  If( CHAR_COUNT .ne. 0 ) Then      ! Char count not 0
! Check for char 1-9
  If( COMMAND .ge. '0' .and. COMMAND .le. '9' ) Then
    DEF_TIME_OUT = Ichar( COMMAND ) - Ichar( '0' )
! Not 1-9 try a command.
  Else
    ERASE_FLAG = .true.      ! Screen erase
    Do 24-1 = 1, COMMAND_MAX
      If( COMMAND .eq. COMMAND_TABLE(1, TABLE_INDEX)) DISPLAY_FLAG = I
24      Continue
    End If
    PURGE = TT_PURGE      !Purge the input buffer next time
  End If

! Get date and time, then dispatch to display routine
25  Call DATE (XDATE)
  Call TIME (XTIME)

  Go to (50,60,99,40) DISPLAY_FLAG

! Refresh the screen (Erase and Redisplay)
30  DISPLAY_FLAG = OLD_DISPLAY      !Redisplay last display
  ERASE_FLAG = .true.
  Go To 25

! Display the HELP frame, set the temporary time-out to 1 minute
40  Type 600, HOME, ERASE_SCREEN      !Display the help frame
  TIME_OUT = 60      !Give the user 1 minute to read it
  DISPLAY_FLAG = OLD_DISPLAY      !When it times out, default old
  ERASE_FLAG = .true.
  Go To 21

! Generate the Status Line for each A/D channel
50  If ( ERASE_FLAG ) Type 300, HOME, ERASE_SCREEN
  Type 100, HOME, XTIME, XDATE
  CHANNEL_COUNT = 0
  Do 51 CHANNEL = 1, MAX_AD_CHANNEL
    If( AD_BLOCK(2, CHANNEL) .ne. 0 ) Then      !If allocated, display info
      Type 200, CHANNEL, STATUS(AD_BLOCK(1, CHANNEL)+1),
      1 (AD_BLOCK(J, CHANNEL), J = 2, 6 )
      CHANNEL_COUNT = CHANNEL_COUNT + 1
    Else      !If not allocated, say so
      Type 900, CHANNEL, '<Unused>', ERASE_LINE
    End If
51  Continue

```

```

PID_COUNT = 0
Do 52 PID_INDEX = 1, MAX_PID
PID = CONNECT_BLOCK(PID_INDEX,1)
If ( PID .ne. 0 ) PID_COUNT = PID_COUNT + 1
52 Continue

Type 400,ERASE_LINE,PID_COUNT,CHANNEL_COUNT
OLD_DISPLAY = DISPLAY_FLAG
ERASE_FLAG = .false.
Go to 20

!
! Status display via process (PID)
!
60 If ( ERASE_FLAG ) Type 300, HOME,ERASE_SCREEN
Type 100, HOME,XTIME,XDATE
PID_COUNT = 0 ! Number of connected processes
CHANNEL_COUNT = 0 ! Number of allocated channels
Do 61 PID_INDEX = 1, MAX_PID
PID = CONNECT_BLOCK(PID_INDEX,1)
If ( PID .ne. 0 ) Then
PID_COUNT = PID_COUNT + 1
OLD_COUNT = CHANNEL_COUNT
Do 62 CHANNEL = 1, MAX_AD_CHANNEL
If( AD_BLOCK( 2,CHANNEL) .eq. PID ) Then !If right PID, display info
Type 200, CHANNEL, STATUS(AD_BLOCK(1,CHANNEL)+1),
1 (AD_BLOCK(J,CHANNEL), J = 2,6 )
CHANNEL_COUNT = CHANNEL_COUNT + 1
End IF
Continue
62 If (OLD_COUNT .eq. CHANNEL_COUNT ) Type 800, '<None>',PID,ERASE_LINE
End IF
61 Continue
Type 400,ERASE_LINE,PID_COUNT,CHANNEL_COUNT,ERASE_SCREEN
OLD_DISPLAY = DISPLAY_FLAG
ERASE_FLAG = .false.
Go to 20

!
! Exit
!
99 Call Exit

!
! Format Statments
!
100 Format(1X,2A1,' Lab IO Status as of ',A,' ',A//
1' Channel Status PID Tics/Sample Buffer Size
1 Buffers '/')

200 Format(15,5x,A8,Z10,4I12)

300 Format(' ',4A1)

400 Format(' ',2A1/' Totals: ',I2,' Processes connected ',I2,' Channels
1 allocated'/')

```



```
500  Format(' '7A1)
600  Format(' '4A1/
1'  The following commands are available://'
1'  VT100  VT52  any'//
1'  -----
1'  PF1    red    C  Channel Display'//
1'  PF2    blue   P  Process Display'//
1'  PF3    grey   H  Help Display'//
1'  PF4    n/a    E  Exit'//
1'  To change display time, type a digit 0-9 for the desired time'//)
700  Format(A)
800  Format(' ',A6,11X,Z10,2A1)
900  Format(15,5x,A8,2A1)
      End
! [End of File]
```


XALINK
MAR

DRMASTER
FOR

XMESSAGE
MAR

LABIOCOM
FOR

LABIOPEAK
FOR

LABIOSTR
COM

XATEST
FOR

LABIDEMO
COM

LABMBXDEF
FOR

LABIOSAMP
FOR

MAILCOMPRESS
COM

LABCHNDEF
FOR

CONNECT
COM

LABIOCON
FOR

LABIDEMO
FOR

PEAK
FOR

DRCOPYBLD
COM

XIDRIVER
MAR

LABIOSEC
FOR

DRSLAVE
FOR

LABIOACQ
FOR

LABIOCOMP
COM

LABIOSTAT
FOR

TESTLABIO
FOR